

# Chapter 7

## High-Throughput Immunofluorescence Microscopy Using Yeast Spheroplast Cell-Based Microarrays

AQI Wei Niu, G. Traver Hart, and Edward M. Marcotte

### Abstract

We have described a protocol for performing high-throughput immunofluorescence microscopy on microarrays of yeast cells. This approach employs immunostaining of spheroplasted yeast cells printed as high-density cell microarrays, followed by imaging using automated microscopy. A yeast spheroplast microarray can contain more than 5,000 printed spots, each containing cells from a given yeast strain, and is thus suitable for genome-wide screens focusing on single cell phenotypes, such as systematic localization or co-localization studies or genetic assays for genes affecting probed targets. We demonstrate the use of yeast spheroplast microarrays to probe microtubule and spindle defects across a collection of yeast strains harboring tetracycline-down-regulatable alleles of essential genes.

**Key words:** Yeast, immunofluorescence, high-throughput microscopy, cell microarrays, microtubule.

---

### 1. Introduction

DNA microarrays and mass spectrometry, among other approaches, have proved to be powerful tools for direct assay of gene function on a genome-wide scale (1, 2). However, such approaches often do not monitor gene and protein behavior in intact cells nor do they assay behavior at the single cell level (3). Nonetheless, it is often desirable to measure single cell phenotypes, and high-throughput microscopy offers a viable strategy for such assays, especially when used in combination with available libraries of genetically distinct cells. To aid in performing such assays, we have described an extension of the spotted cell

Niu, Hart, and Marcotte

49 microarray platform (4) to perform high-throughput immunofluorescence assays on yeast cells. This approach, dubbed yeast  
50 spheroplast microarrays, enables analysis of localization phenotypes of probes of interest (e.g., any target of an affinity reagent  
51 such as an antibody) in thousands of genetically distinct cells in parallel.  
52  
53

54  
55 Yeast spheroplast microarrays have several advantages compared with high-throughput assays in 96- or 384-well plates. First,  
56 the high capacity of the arrays (up to ~5,000 spots per array) minimizes the use of expensive reagents to screen a library of cells by  
57 limiting the use of antibodies or dyes to single microscope slides. Second, the high density of samples in the arrays increases data  
58 acquisition speed by automated microscopy. Third, yeast spheroplast microarrays provide a platform for highly parallel and repro-  
59 ducible large-scale screens. More than a hundred cell microarrays can be manufactured from the same sample sources under the  
60 same conditions, and microarrays can be stored at  $-80^{\circ}\text{C}$  at least for a month without apparent loss of immunostaining signals. As  
61 each slide can in principle be probed with a different antibody, this allows many independent imaging assays to be performed on  
62 cells drawn from the same biological samples.  
63  
64

65  
66 Performing such high-throughput immunofluorescence experiments in yeast requires spheroplasting the cells to allow  
67 antibody access to intracellular contents. One of the primary limitations of high-throughput immunofluorescence imaging  
68 of yeast cells is the difficulty associated with uniformly spheroplasting cells in a strain collection using a single spheroplasting  
69 condition, especially given that optimal treatment times can vary with genetic background. Therefore, spheroplast conditions must  
70 be optimized by testing the time of zymolyase treatment when manufacturing arrays from different yeast cell libraries.  
71  
72

73  
74 To enable such experiments, we present here a protocol for high-throughput spheroplasting of yeast cells and construction  
75 and imaging of spheroplast cell microarrays (*see Fig. 7.1* for an overview of the protocol). As a demonstration of the construction  
76 and utility of such arrays, yeast spheroplast chips were built from the *TetO7* promoter collection covering 75% of essential genes  
77 (5), then probed with anti-alpha tubulin antibody to examine the effects on the organization of microtubule and spindle when  
78 essential genes were downregulated (**Fig. 7.2**).  
79

80  
81 Yeast microtubule organization depends strongly on cell cycle progression (6, 7). In wild-type G1-phase cells, microtubules  
82 form star-like cytoplasmic arrays radiating from the microtubule organizing center (**Fig. 7.3a**). During S/G2 phase, cytoplasmic  
83 astral microtubules reorient towards the bud tip and penetrate the bud (**Fig. 7.3b**). When cells enter mitosis, microtubules  
84 rearrange into a bipolar spindle positioned through the bud neck. In anaphase, cytoplasmic microtubules shorten, the spindle  
85  
86  
87  
88  
89  
90  
91  
92  
93  
94  
95  
96

High-Throughput Immunofluorescence Microscopy




Fig. 7.1. Work flow of the yeast spheroplast microarray protocol. Yeast cells are cultured, fixed, and spheroplasted in 96-well plates, e.g., using an automated liquid handling robot. Subsequently, they are robotically printed onto poly-L-lysine microscope slides using a slotted steel pin-based DNA microarray robot. The resulting slides each contain ~5,000 spots. Each slide can be stored at -80°C or can be probed with a specific antibody immediately after printing. Images are acquired through automated microscopy. The different steps of the protocol are indicated in boxes by arrows. The approximate time required for each step is indicated near the box. Steps at which the protocol can be paused are indicated in the diagram by pairs of short lines across the arrows. Times indicated for each step are approximate and will depend upon the precise equipment used for the experiment.

This figure will be printed in b/w





Fig. 7.2. A microarray printing robot suitable for printing yeast spheroplast microarrays, shown alongside an example microarray (a). A Stanford-style Newport DNA microarray print robot with 12 pins directly delivers yeast spheroplasts from 96-well plate onto slides (a) (4). A yeast spheroplast array carrying ~700 yeast *TetO7*-promoter strains (spot diameter, 200 μm; spot-to-spot distance, 410 μm) is shown in (b).

elongates, and sister chromosomes are segregated (Fig. 7.3c). When cytokinesis is completed, the spindle disassembles (Fig. 7.3d).


In the example spheroplast array, abnormal microtubule structures could be observed in cells with knockdown of genes

Niu, Hart, and Marcotte



145  
146  
147  
148  
149  
150  
151  
152  
153  
154  
155  
156  
157  
158  
159  
160  
161  
162  
163  
164  
165 Fig. 7.3. Examples of microtubule and spindle morphology in different yeast cell cycle stages from (a) to (d) show G1  
166 phase cells, S/G2 phase cells, anaphase cells, and completed cytokinesis, respectively. Nuclei and microtubules are  
167 depicted.

168  
169 whose mutations are known to cause defective spindle forma-  
170 tion. For example, Spc97 and Spc98 are two major compo-  
171 nents of the microtubule nucleating Tub4 complex, involved  
172 in microtubule nucleation and mitotic spindle organization  
173 (8). Downregulation of *SPC97* and *SPC98* resulted in short  
174 spindles, elongated cytoplasmic microtubules, and defective  
175



176  
177  
178  
179  
180  
181  
182  
183  
184  
185  
186  
187  
188  
189  
190 Fig. 7.4. Example images for *TetO7* promoter strains exhibiting defective microtubule and spindle structure (a, b). Down-  
191 regulation of *SPC97* and *SPC98* resulted in short spindles and elongated cytoplasmic microtubules (c). Downregulation  
192 of *TID3* caused defective spindle elongation and uneven nuclear division.

## High-Throughput Immunofluorescence Microscopy

193 chromosome segregation, thus showing a defect resembling that  
 194 of the corresponding temperature-sensitive mutants (Fig. 7.4a,b).  
 195 Similarly, Tid3 is a component of the kinetochore-associated  
 196 Ndc80 complex, involved in chromosome segregation and micro-  
 197 tubule nucleation (9, 10). Downregulation of *TID3* caused  
 198 defects in spindle elongation and nuclear division (Fig. 7.4c).  
 199 Thus, the potential can be seen of spheroplast chips for high-  
 200 throughput immunofluorescence imaging-based genetic screens  
 201 in yeast.

---

## 204 2. Materials

### 207 2.1. Yeast Cell 208 Growth and Fixation

- 209 1. Yeast strains: for the example array shown here, we used  
210 the yeast TetO7-promoter strain collection (5) (Open  
211 Biosystems).
- 212 2. Growth medium: YPED medium (1% yeast extract, 2% pep-  
213 tone, and 2% D-dextrose) supplemented with 200 µg/ml  
214 doxycycline.
- 215 3. 37% Formaldehyde.
- 216 4. 1× PBS buffer: 8 g NaCl, 0.2 g KCl, 1.44 g Na<sub>2</sub>HPO<sub>4</sub>, and  
217 0.24 g KH<sub>2</sub>PO<sub>4</sub> in 800 ml of distilled H<sub>2</sub>O. Adjust pH to  
218 7.4 with HCl. Add H<sub>2</sub>O to 1 l.
- 219 5. Spheroplast buffer: 1.2 M sorbitol, 0.1 M KH<sub>2</sub>PO<sub>4</sub>, pH 7.5.

### 221 2.2. Zymolase 222 Treatment

- 223 1. Zymolase solution: Spheroplast buffer supplemented with  
224 0.025 mg/ml zymolyase 20T (Seikagaku corporation) and  
225 0.4 µl β-mercaptoethanol.

### 228 2.3. Preparation of 229 Poly-L-Lysine Glass 230 Slides

- 231 1. Glass slide wash solution: dissolve 33 g NaOH pellets in  
232 133 ml H<sub>2</sub>O, add additional NaOH pellets until solution  
233 is saturated, and then add 200 ml 95% EtOH.
- 234 2. Poly-L-lysine solution: 250 ml H<sub>2</sub>O, 35 ml 1× PBS (fil-  
235 tered), 50 ml poly-L-lysine.
- 236 3. Glass slide dish.
- 237 4. Glass slide rack.

### 238 2.4. Spheroplast 239 Microarray Printing

- 240 1. DNA microarray printing robot. Here, we use a Stanford-  
style Newport printer.
2. Printer pins. Here, we use conically tapered 1/16 inch diam-  
eter stainless steel printing tips with 0.0015 inch slots (Majer  
Precision Engineering).
3. 0.5× SSC.

Niu, Hart, and Marcotte

241 **2.5. Storage**  
 242 **and Recovery of**  
 243 **Yeast Spheroplast**  
 244 **Microarrays**

1. 95% Ethanol.

245 **2.6. Immunostaining**  
 246 **of Yeast Spheroplast**  
 247 **Microarrays**

1. Methanol,  $-20^{\circ}\text{C}$ .
2. Acetone,  $-20^{\circ}\text{C}$ .
3. Blocking buffer: 3% BSA in  $1\times$  PBS.
4. Primary antibody: Anti-bovine  $\alpha$ -tubulin, mouse monoclonal antibody, and  $4\ \mu\text{g}/\text{ml}$  in blocking buffer.
5. Secondary antibody: Texas-red conjugated goat anti mouse IgG (H+L),  $4\ \mu\text{g}/\text{ml}$  in blocking buffer.
6. VECTASHIELD hard set mounting medium (Vector Laboratories, Inc) with DAPI.
7. Clear nail polish.
8. Staining jars.

248  
 249  
 250  
 251  
 252  
 253  
 254  
 255  
 256  
 257  
 258  
 259 **2.7. Automated**  
 260 **Imaging**

1. Automated fluorescence microscope.

---

261  
 262  
 263  
 264  
 265  
 266 **3. Methods**

267  
 268  
 269 **3.1. Yeast Cell**  
 270 **Growth and Fixation**

1. Grow or treat cells in 96-well plates in  $170\ \mu\text{l}$  proper medium.
2. Fix cells immediately after growth or treatment in proper growth medium with 1/10 volume 37% formaldehyde for 1 h at  $30^{\circ}\text{C}$ .
3. Wash cells with  $1\times$  PBS buffer twice and resuspend cells in  $200\ \mu\text{l}$  spheroplast buffer.
4. At this point, cells can be stored at  $4^{\circ}\text{C}$  overnight or spheroplasted immediately as in [Section 3.2](#).

271  
 272  
 273  
 274  
 275  
 276  
 277  
 278  
 279  
 280 **3.2. Zymolyase**  
 281 **Treatment**

1. Spin down cells and resuspend cells in  $200\ \mu\text{l}$  zymolyase solution (*see Note 1*).
  2. Incubate cells for 2 h at  $30^{\circ}\text{C}$ . This is a critical step and should be monitored closely. *See Fig. 7.5* for images of the desired degree of spheroplasting.
  3. Centrifuge cells at low speeds ( $<3,000\ \text{rpm}$ ) and wash cells twice with spheroplast buffer, and resuspend cells in  $200\ \mu\text{l}$  spheroplast buffer.
- 282  
 283  
 284  
 285  
 286  
 287  
 288

High-Throughput Immunofluorescence Microscopy

289  
290  
291  
292  
293  
294  
295  
296  
297  
298  
299  
300  
301  
302  
303  
304  
305  
306  
307  
308  
309  
310  
311  
312  
313  
314  
315  
316  
317  
318  
319  
320  
321  
322  
323  
324  
325  
326  
327  
328  
329  
330  
331  
332  
333  
334  
335  
336




Fig. 7.5. Examples of morphology changes of yeast cells treated with zymolyase to differing extents. Cells should be a dark, translucent gray after an appropriate degree of digestion (d). (11). Bright cells are insufficiently digested (a–c) (11). “Ghost cells” (pale gray with little any internal structure) have been over-digested (e, f) (11). Spheroplasting protocols for cells typically employ a relatively high concentration of zymolyase and short digestion times (less than 30 min). However, given that yeast spheroplast microarrays are designed for large-scale experiments usually analyzing hundreds or thousands strains at a time, 30 min digestion times are impractical and unforgiving for small errors in timing. We therefore employed lower concentrations of zymolyase and longer digestion times in order to get an appropriate degree of digestion. The optimal concentration for the Tet07 strain collection in our example assay was 0.025 mg/ml; the digestion time was 2 h (d).

**3.3. Preparation of Poly-L-Lysine Glass Slides (Described for 30 Slides)**

4. At this point, spheroplasted cells can be stored at 4°C overnight or printed immediately as in [Section 3.3](#) (see [Note 2](#)).
1. Rinse slide dishes completely with H<sub>2</sub>O in a plastic basin four times to wash away any dust and debris that might be on the slides (see [Note 3](#)).
2. Place racks in glass slide dishes and pour slide wash solution over slides and cover. Shake gently for 2 h (see [Note 4](#)).
3. Rinse each rack with tap-distilled water two times, making sure to rinse both surfaces of each slide.
4. Rinse the racks with double distilled H<sub>2</sub>O at least two times. Make sure that you have removed all of the NaOH.
5. Place one rack of slides into glass slide dishes, draining briefly, and pour poly-L-lysine solution to cover each set of slides. Shake gently for 30 min.
6. Rinse each rack with distilled H<sub>2</sub>O twice in a plastic basin by removing poly-L-lysine from slides.

Niu, Hart, and Marcotte

AQ2 337  
 338  
 339  
 340  
 341  
 342  
 343  
 344  
 345  
 346  
 347  
 348  
 349  
 350  
 351  
 352  
 353  
 354  
 355  
 356  
 357  
 358  
 359  
 360  
 361  
 362  
 363  
 364  
 365  
 366  
 367  
 368  
 369  
 370  
 371  
 372  
 373  
 374  
 375  
 376  
 377  
 378  
 379  
 380  
 381  
 382  
 383  
 384

**3.4. Spheroplast  
 Microarray Printing**

7. Place racks on top of 2 folded over paper towels in a tabletop, centrifuge, and spin at 600 rpm at room temperature for 5 min to dry.
8. Store slides in a clean plastic box at room temperature.

In principle, any slotted pin-based DNA microarray printer can be used. We present specific protocols for the Stanford-style Newport printer, using Arraymaker control software.

1. Depending on the array printer, it may be necessary to set the print depth to be appropriate for 96-well plates prior to loading a plate in the plate holder (*see Note 5*).
2. Use 12 pins, skipping a pin in each direction. Pins are placed in hole numbers 2, 4, 10, 12, 18, 20, 26, 28, 34, 36, 42, and 44 (*see Note 6*).
3. Position poly-L-lysine-coated slides on the printing platter, taping in position if necessary.
4. Fill sonicating water bath with 0.5×SSC.
5. For plate alignment, due to the variability in well volumes, align printer pins so that the liquid level is just below the start of the beveled portion of the printing tip.
6. Align the slide position as typical for DNA microarray printing.
7. Set wash time at 6,000 ms, dry time at 8,000 ms, load time at 1,000 ms, the spot-to-spot distance as 410 μm. Use “slow pickup.”
8. Test print by using a test plate with the wild-type samples before the large-scale printing, and check the cell counts under a microscope (*see Note 7*).
9. Start the formal print runs once the test print results are satisfying.
10. When the print is complete, tape slides into a swinging bucket plate centrifuge (max: 5 slides/bucket) or use custom-made plastic flat slide holders, and spin flat at 1,500–2,000 rpm for 5 min.

**3.5. Storage  
 and Recovery of  
 Yeast Spheroplast  
 Microarrays**

1. Store spheroplast chips at 4°C for overnight or at –80°C for long-term storage (**Fig. 7.6**) (*see Note 8*).
2. To recover a slide after storage at –80°C, fill a 50 ml Falcon tube with 95% EtOH. Dunk slide into the tube immediately upon taking slide out of the freezer (*see Note 9*).
3. Take the slide out quickly and put it into a clean and dry 50-ml Falcon tube, and then dry the slide by spinning in an empty 50-ml tube without cap at 600–700 rpm for 5 min.

High-Throughput Immunofluorescence Microscopy

385  
386  
387  
388  
389  
390  
391  
392  
393  
394  
395  
396  
397  
398  
399  
400  
401  
402  
403  
404  
405  
406  
407  
408  
409  
410  
411  
412  
413  
414  
415  
416  
417  
418  
419  
420  
421  
422  
423  
424  
425  
426  
427  
428  
429  
430  
431  
432




Fig. 7.6. Spheroplast microarrays can be stored at least a month at  $-80^{\circ}\text{C}$  without substantial effects on cellular morphology or subcellular structure, as shown here for microtubule structure. Compared to fresh chips (a), chips stored at  $-80^{\circ}\text{C}$  still maintain intact microtubule structures after one month (b). Nuclei and microtubules are depicted.

**3.6. Immunostaining of Yeast Spheroplast Microarrays**

1. Permeabilize cells in cold methanol for 6 min and then in cold acetone for 30 s. Air-dry slides completely.
2. Wash slides for 5 min with  $1\times$  PBS.
3. Place  $\sim 200\ \mu\text{l}$  blocking buffer on the slide. Trim a piece of parafilm to the size of the slide and gently lay over the sample. Incubate for 30 min at  $30^{\circ}\text{C}$  in a sealed humidity chamber (*see Note 10*).
4. Carefully remove the parafilm and add  $\sim 200\ \mu\text{l}$  primary antibody onto the slide. For the anti-bovine  $\alpha$ -tubulin mouse monoclonal antibody used as an example (Invitrogen, cat. no. A-11126). Incubate for 2 h at  $30^{\circ}\text{C}$ .
5. Wash the slide three times in staining jars with  $1\times$  PBS, 5 min each time.
6. Add  $\sim 200\ \mu\text{l}$  secondary antibody onto the slides. Incubate for 2 h at  $30^{\circ}\text{C}$  and in the dark (*see Note 11*).
7. Wash the slide three times in staining jars with  $1\times$  PBS, 5 min each time.
8. Add  $\sim 60\ \mu\text{l}$  VECTASHIELD hard-set mounting medium and cover with cover slip.
9. Let the mounting medium dry out for 15 min at room temperature or overnight at  $4^{\circ}\text{C}$ . Seal the slides with clear nail polish.
10. Store the slides at  $4^{\circ}\text{C}$  and protect the slides from light. Immunostained yeast spheroplast microarrays can be stored at  $4^{\circ}\text{C}$  for at least a month with no apparent loss of signal.

AQ3

Niu, Hart, and Marcotte

### 3.7. Automated Imaging

In principle, any automated fluorescent microscope can be used to capture cell chip images. Here, we have presented a protocol specific to the Nikon TE2000 inverted epifluorescence microscope with NIS Elements controller software.

1. Turn on the microscope and all components. Initialize the stage. With a dummy slide, calibrate the stage and slide mounting bracket to ensure a flat focal plane across the sample area on the slide. This can be accomplished by marking one side of the slide with reference marks in the corners of the slide (using, for example, a fine-point Sharpie marker), navigating to each mark, and adjusting the calibration screws in the mounting plate until all reference spots can be viewed clearly at the desired magnification by translating X and Y without adjusting focus.
2. Edit the `create_block_coordinates` macro script (*see* **Note 12**) to reflect the block size on the printed microarray. Block size is the number of rows and columns of spots printed by a single pin on the array. Also specify the spot pitch—the distance between the printed spots in the array construction (in this protocol, 410  $\mu\text{m}$  was used in lines 24 and 25 of the script).
3. Mount the slide and manually scan to the first printed spot in the first block. Focus on the spot in DIC at the desired magnification. 60 $\times$  or greater magnification is recommended.
4. Execute the `create_block_coordinates` macro. An ND Experiment file will be created that specifies the coordinates of all spots in the current block.
5. In the ND Experiment dialog box, specify the desired immunofluorescence wavelengths to be captured. Also specify whether to autofocus on every spot. This should be unnecessary if the stage and slide mounting bracket are well calibrated (*see* **Section 3.7**, Step 1). Finally, specify the output file.
6. Execute the ND Experiment. The microscope will image the current block.
7. Repeat Steps 3–6 for each block. A 12-pin print will have 12 blocks.

AQ4

---

## 4. Notes

1. The spheroplast buffer with zymolyase and  $\beta$ -mercaptoethanol should be prepared fresh before use.

### High-Throughput Immunofluorescence Microscopy

481  
482  
483  
484  
485  
486  
487  
488  
489  
490  
491  
492  
493  
494  
495  
496  
497  
498  
499  
500  
501  
502  
503  
504  
505  
506  
507  
508  
509  
510  
511  
512  
513  
514  
515  
516  
517  
518  
519  
520  
521  
522  
523  
524  
525  
526  
527  
528

2. Note that yeast cells without cell walls are quite fragile, and in general spheroplasted cells should be handled carefully and centrifuged at low speeds.
3. We recommend wearing powder-free gloves at all times, as talcum powder will stick to slides and leave visible debris in the resulting images.
4. This wash removes residual oil and debris from slides. The wash solution should be prepared fresh before use.
5. On Stanford-style printers optimized for DNA microarray printing, this step is important, because 96-well plates are often taller than 384-well plates. If the vertical printing depth (Z) is not set high enough, the printing pins will hit the 96-well plate bottom.
6. If using the Stanford-style arrayer and Arraymaker software, be sure to use the 48-tip printing designation in the Arraymaker software.
7. Ideally, there should be ~20–30 cells per spot. If cell count is low, consider resuspending cells in smaller volumes.
8. Empirically, we find that spheroplast arrays are fine if stored at 4°C for a few days, but for longer periods, it is better to store them at –80°C. Storage times up to 1 month have been tested.
9. This step is to prevent condensation on the slide. Condensation will smear the spots, cause potential cross-contamination, and make it very difficult to grade the spots on the slide.
10. It is important to not let the slides dry out from this step on.
11. Keep slides in the dark from this point on to minimize photobleaching of fluor.
12. The create\_block\_coordinates script referenced in protocol [Section 3.7](#), Step 2-on, suitable for use with NIS Elements microscope controller software (specifically tested on NIS Elements AR (Advanced Research) v3.0.)

AQ5

```
//Automatic Spot Coordinate calculation for a spotted
cell microarray
//Program moves in NORMAL PRINT PATTERN
//User inputs x,y,z at first, last spot
//Program calculates x,y,z of every spot
//by Traver Hart, 1/22/07

// BE SURE STAGE IS INITIALIZED BEFORE RUNNING THIS
MACRO
char Debug[1024];
```

Niu, Hart, and Marcotte

```

529     int main()
530     {
531
532         double dYSpotBegin,dXSpotBegin, dZSpotBegin;
533         double dYSpotEnd,dXSpotEnd, dZSpotTwo, dZSpotThree;
534         double dYOffset,dXOffset,dyTemp,dxTemp, dzOffset,
535             dZYOffset, dZXOffset, dzTemp;
536         int i,j;
537         int numRows, sectorWidth;
538         sectorWidth = 21; // (y direction)
539         numRows = 21;     // (x direction)
540
541         //Set distance between spots - specified as 410
542         microns in arrayer
543         dXOffset=410;
544         dYOffset=410;
545
546         Int_Question("Automated Print Block ND
547         Initialization","Please move to the TOP
548         LEFT spot (should be position 1,1 of block).","OK"
549         , "", "", "", 1,1);
550         StgGetPos (&dXSpotBegin,&dYSpotBegin,&dZSpotBegin);
551         dzTemp = dZSpotBegin;
552
553         for(i=0;i<numRows;++i)
554         {
555             for(j=0;j<sectorWidth;++j)
556             {
557                 dxTemp=(dXSpotBegin - (dXOffset*i) );
558                 dyTemp=(dYSpotBegin - (dYOffset*j) );
559                 ND_AppendMultipointPoint(dxTemp,dyTemp,
560                 dzTemp,"sample");
561             }
562         }
563         _ND_CreateAndRunExperiment();
564     }
565
566
567
568
569
570
571

```

---

## Acknowledgments

This work was supported by grants from the N.S.F., N.I.H., and Welch (F-1515) and Packard Foundations to E.M.M.

## High-Throughput Immunofluorescence Microscopy

## References

577  
578  
579  
580  
581  
582  
583  
584  
585  
586  
587  
588  
589  
590  
591  
592  
593  
594  
595  
596  
597  
598  
599  
600  
601  
602  
603  
604  
605  
606  
607  
608  
609  
610  
611  
612  
613  
614  
615  
616  
617  
618  
619  
620  
621  
622  
623  
624

1. Lockhart, D. J., Winzeler, E. A. (2000) Genomics, gene expression and DNA arrays. *Nature* 405, 827–836.
2. Aebersold, R., Mann, M. (2003) Mass spectrometry-based proteomics. *Nature* 422, 198–207.
3. Levisky, J. M., Shenoy, S. M., Pezo, R. C., Singer, R. H. (2002) Single-cell gene expression profiling. *Science* 297, 836–840.
4. Narayanaswamy, R., Niu, W., Scouras, A. D., Hart, G. T., Davies, J., Ellington, A. D., Iyer, V. R., Marcotte, E. M. (2006) Systematic profiling of cellular phenotypes with spotted cell microarrays reveals mating-pheromone response genes. *Genome Biol* 7, R6.
5. Davierwala, A. P., Haynes, J., Li, Z., Brost, R. L., Robinson, M. D., Yu, L., Mnaimneh, S., Ding, H., Zhu, H., Chen, Y., et al. (2005) The synthetic genetic interaction spectrum of essential genes. *Nat Genet* 37, 1147–1152.
6. Winsor, B., Schiebel, E. (1997) Review: an overview of the *Saccharomyces cerevisiae* microtubule and microfilament cytoskeleton. *Yeast* 13, 399–434.
7. Kilmartin, J. V., Adams, A. E. (1984) Structural rearrangements of tubulin and actin during the cell cycle of the yeast *saccharomyces*. *J Cell Biol* 98, 922–933.
8. Knop, M., Pereira, G., Schiebel, E. (1999) Microtubule organization by the budding yeast spindle pole body. *Biol Cell* 91, 291–304.
9. Zheng, L., Chen, Y., Lee, W. H. (1999) Hec1p, an evolutionarily conserved coiled-coil protein, modulates chromosome segregation through interaction with SMC proteins. *Mol Cell Biol* 19, 5417–5428.
10. Wigge, P. A., Kilmartin, J. V. (2001) The *ndc80p* complex from *Saccharomyces cerevisiae* contains conserved centromere components and has a function in chromosome segregation. *J Cell Biol* 152, 349–360.
11. Burke, D., Dawson, D., Stearns, T. (2000) *Methods in Yeast Genetics*. Cold Spring Harbor Laboratory Press, Cold Spring Harbor, NY.

625  
626  
627  
628  
629  
630  
631  
632  
633  
634  
635  
636  
637  
638  
639  
640  
641  
642  
643  
644  
645  
646  
647  
648  
649  
650  
651  
652  
653  
654  
655  
656  
657  
658  
659  
660  
661  
662  
663  
664  
665  
666  
667  
668  
669  
670  
671  
672

## Chapter 7

---

Q. No.	Query
AQ1	Please provide e-mail id for the contributors “Wei Niu and G. Traver Hart”.
AQ2	Please check sentence “Place racks on top...” for clarity.
AQ3	Please check sentence “For the anti-bovine $\alpha$ -tubulin. . . . .” for clarity.
AQ4	“Ster 3.7.1” has been changed to “Section 3.7, Step 1,”. Please check.
AQ5	“Ster 3.7.2” has been changed to “Section 3.7, Step 2.” Please check.

---

UNCORRECTED PROOF